
Subject: Infection Control, COVID-19 and Support Surface Cleaning.

COVID-19 and its quick spread has brought infection control to the forefront of community and facility healthcare. This has prompted CMS to direct surveyors to place more focus on infection control policies and procedures when conducting site visits in post-acute care facilities. ¹

Infection control has always been a critical—though sometimes not standardized—component of healthcare. This is especially true in post-acute care settings, where smaller staff, more staff turnover and training challenges can result in reduced adherence to infection control policies and procedures. ²

Post-acute care providers are also challenged by a patient population that in general is more debilitated than in acute care. These patients tend to have more illnesses, co-morbidities and symptoms, which make it more difficult to recognize small changes in their condition. Additionally, many post-acute patients cannot effectively communicate with staff when they have new health issues.

Considered together, these multiple challenges help explain why post-acute facilities have been described as “Ground Zero” in the battle against the spread of COVID-19.

PressureGuard® Support Surfaces: Designed with infection control in mind.

By their nature, typical “roll-up”-style air therapy surfaces found in many facilities can make effective in-facility cleaning a challenge. In contrast, the stable outer shell design of Span-America’s air therapy surfaces allows them to be cleaned and disinfected in place much like a standard replacement mattress. This enables a level of thorough, effective in-room cleaning and disinfection that is extremely difficult to achieve with roll-up surfaces.


In addition, Span uses polycarbonate-fortified *LifeSpan*® premium stretch fabrics on several of its air surfaces. The polycarbonate additive in the coating of *LifeSpan* fabric provides a level of durability in the presence of strong chemical cleaners and disinfectants unachievable with the standard stretch fabrics used on most typical surfaces. This protection also reduces the potential for fluid ingress, which often results from chemical degradation of a fabric’s fluid-resistant coating over time. ³

In support of facility infection control efforts, Span-America recommends the following:


- Use an approved antimicrobial agent.
 - [For all the registered antimicrobial agents for disinfecting, please refer to List N: EPA’s Registered Antimicrobial Products for Use Against Novel Coronavirus SARS-CoV-2, the Cause of COVID-19.](#)
- Wipe down in place and adhere to the suggested manufacturer’s wait/dwell time.
- Properly rinse all cleaning agents and disinfecting chemicals after application with clean water and allow to dry.

- Over time, even approved cleaning solutions may cause damage to the integrity of the fabrics used on medical mattresses. Cleaning agents that are strong enough to be efficacious cleansers and disinfectants may cause degradation of the same fabrics on which they are being used. To minimize negative impact and maximize the usable life of the medical mattress covers, attention needs to be paid to the dwell time, or time the when the cleaning solution is in contact with the fabric. This must be monitored and kept to the required time specified on the manufacturer's instructions.
- For PressureGuard Surfaces that include microclimate management, the outer coverlet can be removed and machine laundered. Air dry only.
- During the cleaning process, covers should be thoroughly inspected for wear, cuts and tears, or damage of any kind. If the cover has been compromised in any way, the product should be removed from service and cover replaced. If the internal foam/system appears compromised, i.e., discolored, damaged, wet, or has odor, the entire product should be disposed of immediately.
- Terminal cleaning protocol between patients should be strictly followed.

Complete cleaning instructions can be found at spanamerica.com.



Kathy Whittington, RN, MSN, CWCN
Manager, Clinical Services



Derrill Darby,
Director, Quality

1. <https://www.cms.gov/newsroom/press-releases/cms-announces-actions-address-spread-coronavirus>
2. <https://www.mcknights.com/news/workforce-crisis-among-reasons-why-nursing-homes-are-most-at-risk-for-covid-19>
3. Span-America.com: [LifeSpan Fact sheet](#).

For more information:

<https://www.cdc.gov/hai/index.html>

https://www.epa.gov/sites/production/files/2020-03/documents/sars-cov-2-list_03-03-2020.pdf

PS Infection Control, COVID-19 and Support Surface Cleaning
CO# 2334