


ADVANTAGE

— WEIGH SCALE BED —

Span's Weigh Scale low-bed design with advanced clinical positioning and ReadyWide® integrated expansion now available with built-in scale



Position shown:
Auto Contour

The Advantage ReadyWide Weigh Scale Bed with advanced features including mobility at any height, easy optional sleep surface expansion to 42", a multi-function carrier system and a safe working load (SWL) of 500 lbs (227 kg).

Now available with choice of two weigh scale package options:



Dual Controllers: Tethered staff control / lockout pad and separate handheld weigh scale pendant.

OR



Premium Digital Touchscreen with fully integrated staff controls, lockout and weigh scale functions.

