

# PressureGuard® APM<sup>2</sup> Safety Supreme



## Therapy Provided

- ✓ Powered Pressure Redistribution
- ✓ Alternating Pressure
- ✓ Lateral Rotation
- ✓ Raised Safety Perimeter

Combines all the features of the top-selling PressureGuard® APM<sup>2</sup>, plus the added safety benefits of raised perimeter edges. The Safety Supreme incorporates the proprietary Geo-Mattress with Wings® raised perimeter design to help prompt the user away from the edge of the bed. Provides protection from pressure injuries, falls, fire, restraints and siderail entrapment – all in one easy-to-use package.

## QUICK FACTS

### Indications:

(in addition to APM<sup>2</sup> p. 8):

- Patients at risk for falls

### Exclusive Raised Perimeter

*incorporates the unique "edge-to-edge" design to help prompt user away from edge of bed without compromising skin integrity.*



Item #	Description
SAF5875-29	Safety Supreme, 75"L x 35"W x 7"H center, 9"H side
SAF5880-29	Safety Supreme, 80"L x 35"W x 7"H center, 9"H side
SAF5884-29	Safety Supreme, 84"L x 35"W x 7"H center, 9"H side
* 37890	Safety Supreme, 80"L x 42"W x 7"H center, 7"H side
* 33530	Safety Supreme, 84"L x 42"W x 7"H center, 7"H side
5900	Digital control unit

\* Allow up to 10 working days lead time for this size.

**Medicare Approved:** E0277 (Group II)

**Weight limit:** 350 lbs.

**Warranty:** 18 months (all components)

**Mattress Weight:** approx. 20 lbs.

**UL listed:** Yes

All models are certified for conformance to NFPA 101 (Life Safety Code) – ASTM E1590; Cal Tech Bulletin #117 and 16 CFR Parts 1632 and 1633 by an independent testing organization.

