


The "Soft-Pedal"[™]
Low-Profile Cushion
 With *SynerGel[™]

Just 1 7/8" high, this cushion is the ideal choice for wheelchair users who require – or prefer – a low-profile design. Facilitates self-propelled users who wish to have full foot contact with the floor, especially in hemi-height chairs.

The Soft-Pedal[™] cushion's multi-layer design provides comfort and support with the following features:

- a soft memory-foam top
- silicone gel (SynerGel[™]) pad insert for added support and protection
- high-density foam base
- wipe-clean cover with non-slip top and non-skid bottom
- 250 lb. weight limit
- 18-month warranty



**SynerGel is a registered trademark of Synergy Rehab Technologies, Inc.*

Sizes available	Part number◆
16" wide x 16" long	86344 (6/cs)
18" wide x 16" long	82658 (6/cs)
20" wide x 18" long	89156 (6/cs)

◆ add "-1" suffix for "each" quantity