




# Geo-Mattress® UltraMax™ for Stretchers

**Category:**  
Premium Pressure Redistribution  
Support Surface




## Design Advantages Inside...

Engineered, multi-layer  
support core

3-Dimensional Zoning™ design provides anatomically-calibrated levels of support and immersion throughout.



Extra-soft heel section facilitates immersion to provide added heel protection.



UHP 500™ Ultra High Performance foam in high-fatigue seat section; 55% higher support factor yields superior resiliency and recovery.



Geo-Matt® non-shearing surface creates individually articulating cells; Ring-of-Air® channels dissipate heat and moisture.

Now made with



bleach-resistant fabric cover  
(see page 67 for details)

### Specifications

[All models listed are 5" high to facilitate maximum immersion]

| Item #  | Dimensions, overall (L)ength and (W)idth | (H)ead-end taper length | (F)oot-end taper length |
|---------|--|-------------------------|-------------------------|
| 47318   | 75"L x 26"W                              | 7"                      | 0"                      |
| 49297   | 75"L x 26"W                              | 13"                     | 4"                      |
| 97107   | 76"L x 26"W                              | 10"                     | 4"                      |
| 10635   | 76"L x 26"W                              | 12"                     | 5"                      |
| 2879529 | 76"L x 27"W                              | 11"                     | 5"                      |
| 46234   | 76"L x 29"W                              | 10"                     | 8"                      |
| 79012   | 76"L x 30"W                              | 13"                     | 4"                      |
| 98969   | 78"L x 27"W                              | 12"                     | 12"                     |
| 81921   | 2-pc "Gynie"                             | n/a                     | n/a                     |

Other dimensions also available upon request.


**Supports:** up to 700 lbs.

**Warranty:** 3-year, non-prorated

**Mattress Weight:** approx. 20 lbs.

**All components:** 100% latex-free and 100% PBDE-free

**Flammability:** All models are certified for conformance to NFPA 101 (Life Safety Code) – ASTM E1590; Cal Tech Bulletin #117 and 16 CFR Parts 1632 and 1633 by an independent testing organization.



"Universal" hook & loop pattern. Large areas of both hook and loop on bottom of mattress help in anchoring it to a wide variety of stretcher models.

