
## **Encore Structural Highlights**



## **Optional Accessories:**

- Bed extenders to 84"
- Ergonomic 3-position assist device
- Bed & deck angle indicator
- Snap-in pendant holder fits all rails
- Molded plastic half-rail
- Headboard mount pendant holder
- Transport dollies
- Portable backup battery
- Regular backup battery
- · Drainage hooks
- IV pole
- Trapeze

| Height Floor to Mattress Support    | 7-7/8" - 31"              |  |
|-------------------------------------|---------------------------|--|
| Deck Angles - Back Elevation        | Back Elevation 0 – 70°    |  |
| Deck Angles - Knee Elevation        | 0 – 30°                   |  |
| Bed Length - Overall Length         | 82", 86", 90"             |  |
| Bed Length - Retraction Travel      | 12"                       |  |
| Bed Length - Sleep Surface          | 76", 80", 84"             |  |
| Bed Width - With Side Rails         | 40-1/2", 43-1/2", 46-1/2" |  |
| Bed Width - Sleep Surface           | 35", 39", 42"             |  |
| Weight of Bed (without accessories) | 253 lbs / 115 kg          |  |
| Safe Working Load                   | 600 lbs / 272 kg          |  |

| Item #   | Description  |
|----------|--|
| MC9RW    | Encore bed with GlideAlign retractability, with readywide expanders up to 42" wide                             |
| MC9RWZ   | Encore bed with GlideAlign retractability, with readywide expanders up to 42" wide and built-in battery backup |
| R9938-L  | Three position pivot assist - left   |
| R9938-R  | Three position pivot assist - right  |
| R9939    | Three position pivot assist - set  |
| R9308H   | Half head end side assist/rail - set   |
| R9308H-L | Half head end side assist/rail - left  |
| R9308H-R | Half head end side assist/rail - right   |
| Q6690    | Extender kit 4" - to 84"   |
| Q6699    | Dolly set  |
| Q6730    | Universal Trapeze bar  |
| Q6731    | Encore Trapeze mounting frame  |
| Q6708    | Roller bumper kit  |
| Q6703    | Deck angle indicators with brackets  |



For ordering and technical support on beds, accessories and furniture, call: 800-268-8671