


Advanced Clinical Positioning



Trendelenburg



Reverse Trendelenburg



Trapezoidal slat profile provides unmatched strength, ventilation, and ease of cleaning.

Position shown: Comfort Chair

Standard hand pendant

Comes standard with underbed lighting

4" bed extenders to 84" (2134 mm)


Choice of 6 finishes in 8 profiles



Standard Melamine Profiles



Elite HPL Profiles



Height Floor to Mattress Support	8-5/8" - 32-1/4"
Deck Angles - Back Elevation	0 - 70°
Deck Angles - Knee Elevation	0 - 30°
Bed Length - Overall Length	82", 86", 90"
Bed Length - Sleep Surface	76", 80", 84"
Bed Width - Sleep Surface	35", 39", 42"
Weight of Bed (without accessories)	297 lbs / 135 kg
Safe Working Load	500 lbs / 227 kg
Warranty	Mechanical: 5 years Frame: 10 years Welds: 15 years

Item #	Description
MC7WSTML	Weigh Scale Advantage ReadyWide with mobility, floor lock, premium foot end controls and touchscreen
MC7WSML	Weigh Scale Advantage ReadyWide with mobility, floor lock and foot end controls
R9938-R	3-position molded pivot assist rail with right bracket
R9938-L	3-position molded pivot assist rail with left bracket
R9939	Set of 3-position molded pivot assist rails
R9308H	Set of molded head end half side rails
R9308H-L	Left side molded head end half side rail
R9308H-R	Right side molded head end half side rail
R9519	Split side rail, 3/4 bed coverage
Q6770	Battery back up - factory installed
Q6740	Floor replacement arm
Q6741	Retractable power cord cable reel and placement arm
Q6690	Bed extender to 84"
Q6730	Patient helper, trapeze bar, chain and hand grip
Q6731	Trapeze bed mounting bracket
Q6703	Bed and deck angle indicator
Q6699	Bed transportation system
Q6531	Honeycomb roller bumper kit

