

ADVANTAGE

READYWIDE SERIES

Span's low-bed design with advanced clinical positioning and

Readywide® Integrated Expansion



In addition to the hand pendant controller, a foot-end control with staff-only lockout is standard.



Tool-less, "two-click" expansion



Optional moulded quick release plastic half-head assist rails.



Optional, moulded quick release plastic half-head assist rails.

Standard features include:

- Durable, anti-microbial polyester powder-coating of metal components
- Smooth deck design for easier cleaning
- Mobility at any height
- Foot end controller with lockout
- Vascular foot rest with easy ratchet operation
- Beds move vertically with virtually no horizontal motion
- Pendant bed controller

Optional features:

- Advanced mobility as well as floor level locking
- Bed carrier for streamlined bed logistics and maintenance
- Stylish head and foot board options
- Bed extenders

Height Floor to Mattress Support	7-7/8" - 31"
Deck Angles- Back Elevation	0 - 70°
Deck Angles - Knee Elevation	0 - 30°
Bed Length - Overall Length	82", 86", 90"
Bed Length - Sleep Surface	76", 80", 84"
Bed Width - Sleep Surface	35", 39", 42"
Weight of Bed (without accessories)	259 lbs / 118 kg
Safe Working Load	500 lbs / 227 kg 650 lbs / 295 kg
Warranty	Mechanical: 5 years Frame: 15 years Welds: 15 years

Also Available

NEW!

"PLUS" MODEL

650 LBS
SAFE WORKING LOAD

Advanced Clinical Positioning



Item #	Description
MC7W	Advantage ReadyWide with foot end controls
MC7WM	Advantage ReadyWide with mobility & foot end controls
MC7WL	Advantage ReadyWide with floor lock & foot end controls
MC7WML	Advantage ReadyWide with mobility, floor lock & foot end controls
MC7WPML	Advantage ReadyWide Plus (650 lbs SWL) with mobility, floor lock & foot end controls
R9938	3 position molded pivot assist rail
R9939	Set of 3 position molded pivot assist rails
R9308H	Set of molded head end half side rails
R9308H-L	Left side molded head end half side rail
R9308H-R	Right side molded head end half side rail
R9518	Split side rail, 3/4 bed coverage
Q6749	Under bed lighting - factory installed
Q6701	Battery back up - factory installed
Q6770	Battery back up (Plus model) - factory installed
Q6740	Floor placement arm
Q6741	Retractable power cord cable reel & placement arm
Q6690	Bed extender to 84"
Q6730	Patient helper, trapeze bar, chain and hand grip
Q6731	Trapeze bed mounting bracket
Q6703	Bed and deck angle indicator
Q6699	Bed transportation system
Q6531	Honeycomb roller bumper kit



For ordering and technical support on beds, accessories and furniture, call: **800-268-8671**