

## **Span Seating Products**

No wound care protocol or patient care plan is complete until seating needs are addressed. So, whether the goal is complete pressure removal at the injury site, chronic pain relief, therapeutic positioning, pressure injury prevention, or simple comfort, Span America offers a variety of innovative seating solutions to achieve it.



Clinically-proven Geo-Matt®

segmented top surface provides pressure management and shearing control. It also aids in the dispersion of heat and moisture through exclusive Ring-of-Air® channels at the base of each independently articulating segment.

**Featured on:** Geo-Matt Contour, Geo-Matt Wedge, Geo-Matt PRT, Geo-Matt 4" cushion.

Positioning features -

including "bucketing" to
hold the pelvis in place
and contoured channels
to maintain the thighs in
parallel alignment—help the user
maintain a more healthful upright posture
while playing an important role in fall prevention.

**Featured on:** Geo-Matt Contour, EZ-Dish, Equalizer, Isch-Dish, Isch-Dish Thin.

## **Elimination** at the injury site is a key to injury healing. Since

**Pressure** 

injury healing. Since 1989, the patented, pressure-eliminating pocket design of the



Isch-Dish has provided the only ready-made solution for completely offloading sitting pressure from the ischials and coccyx.

Featured on: Isch-Dish

## **Durafilm™ Coating**

a proprietary method of fusing a durable, waterproof barrier to foam increases the longevity of the product.

**Featured on:** Isch-Dish, Isch-Dish Thin.

45



RV#19E