

Seating

Category: Basic Prevention and Comfort

Soft-Pedal™ Low-Profile Cushion

Multi-layer design with a memory-foam top, silicone gel (SynerGel™) pad insert and high-density foam base, provides comfort and support, and just 1-7/8" high. Ideal for self-propelled wheelchair users, or those who require or prefer a low-profile cushion. Includes a wipe-clean cover with a non-slip top and non-skid bottom.




Item #	Description	Qty/Cs
86344	16"Wx16"L, w/cover	6
86344-1	16"Wx16"L, w/cover	1
82658	18"Wx16"L, w/cover	6
82658-1	18"Wx16"L, w/cover	1
89156	20"Wx18"L, w/cover	6
89156-1	20"Wx18"L, w/cover	1

QUICK FACTS RISK LEVEL
 Low Moderate High

Use for:
 • Comfort

Cross Cut Cushion

Economical prevention and comfort. Provides pressure dispersion and patient comfort through segmented top surface. Cover consists of anti-slip top surface and non-skid bottom. Built-in sleeve in bottom of cover allows use of molded Sling-Fill base (not included-see below); 1-year warranty. Height 3".



Item #	Description	Qty/Cs
CC1616-96	Cross-cut Cushion w/ anti-slip cover, 16"W x16"L	6
CC1816-96	Cross-cut Cushion w/ anti-slip cover, 18"W x16"L	6
CC1818-96	Cross-cut Cushion w/ anti-slip cover, 18"W x18"L	6
CC2018-96	Cross-cut Cushion w/ anti-slip cover, 20"W x18"L	6

Fire barrier versions can be ordered to conform with 2000 NFPA 101 (Life Safety Code) and Cal TB# 133. To order, replace the "-96" suffix with "-29" in item number.

QUICK FACTS RISK LEVEL
 Low Moderate High

Use for:
 • Comfort

Convoluted Comfort Cushions

Convoluted comfort cushions are made from high-quality polyurethane foam. Specifications vary to meet differing levels of discomfort.

Note: Convoluted cushions are effective as a comfort measure only, and are not intended for use by patients at risk for pressure ulcers.

Item #	Description	Qty/Cs
50746-010	18" x 16" x 4"P (.75" base) Premium Density	12
50077-109	16" x 16" X 3"P (.5" base) Standard Density	12
50077-112	16" x 16" x 3"P (1.5" base) Standard Density	12
50746-025	18" x 16" x 4"P (.75" base) Standard Density	12
SP663-000	18" x 16" x 4"P (1.5" base) Standard Density	12

Category: Specialty Seating Accessories

Sacral Dish™

Pressure elimination, positioning, and increased comfort for those with, or at high risk, for sacral and coccygeal pressure injuries. Provides pressure-free, ventilated protection at sacrum. Semi-rigid back base and unique shape create supportive foundation to fill in chair's back upholstery. Designed to work in concert with Isch-Dish cushions. 2 year warranty.



Item #	Description	Qty/Cs
321R	Regular 16"Wide x 16"Tall	1
321W	Wide 18"Wide x 16"Tall	1
321RT	Regular Tall 16"Wide x 22"Tall	1
321WT	Wide Tall 18"Wide x 22"Tall	1

QUICK FACTS

Primary Use:

- Sacrococcygeal injury healing or prevention in sitting

Key Advantages:

- Totally off-loads and ventilates sacral area
- Can be used with any cushion
- Especially effective in reclined chair position

Sling-Fill™ Base

Molded ABS plastic insert fills in the gap created by fabric slings. Creates a flat, rigid foundation where needed. One size fits all. Included standard with all Isch-Dish cushions. Fits into sleeve sewn into all covers on cushions with "-96" suffix.



Item #	Description	Qty/Cs
2SFB	Sling-Fill Base, single	1
2SFBN-06	Sling-Fill Base, set of 6	6

(SynerGel is a registered trademark of Synergy Rehab Technologies, Inc.)

