

Span+Aids® Patient Positioners

Category: Extremities



Foot Drop Stop - Standard

Elevates heels, supports ankle position.
Adjustable strap can accommodate ankle contracture, protect against foot drop.
Square base provides stability. Ring-of-Air® design provides ventilation.
50638-511 *6/cs



Limb Elevator

Elevates casted or uncasted leg. Can be used in both supine and side lying positions. 50632-508 *6/cs

50632-308 °6/CS



Foot Drop Stop – Original

Similar to standard version (above), except hexagonal base maximizes comfort by allowing user to lie with leg in multiple positions.


55555-01 *6/cs 50638-311 6/cs



Flared Leg Extension

Elevates leg and knee in extended position.


SP291-000 2/cs



Cradle Boot – Original

Floats heel, supporting foot in neutral position, prevents external leg rotation.

55555-02 *6/cs



Leg Trough

Immobilizes leg and knee in neutral position. Base is contoured to posterior leg.

SP033-000 2/cs



Cradle Boot - Standard

Floats heel, supporting foot in the neutral position. Well-ventilated boot is easy to use and apply. Square base provides stability, discouraging leg rotation. 50644-515 *6/cs



Upper Limb Positioner

Positions the entire upper arm with sloped groove for stability and comfort. Initially designed for use with burn patients.

SP121-000 2/cs



Cock-Up Splint

Designed to help maintain hand and wrist in a desired position. Holds fingers and thumb opposed to prevent contractures. Hand stirrup is adjustable. 50640-512 *8/cs



*Indicates items vacuum packed to save storage. All items conform to California Technical bulletin 117.

52 RV#19E