

Portable Ceiling Lifts

Portable lifts for Mobility at Home

- Designed for easy, intuitive use.
- Lithium-ion battery for optimal performance, less charging time and extended life (up to five years).
- SureClip sling clip for added safety.
- Paired with a gantry system for simple lifting from wheelchair to bed.

For more details about Savaria Ceiling Lifts, specifications, accessories and where to purchase, visit: ceiling-lift.com

SPECIFICATIONS

Product weight	Portable: 11 lb Fixed model: 17 lb Fixed motorized model: 19 lb
Safe Working Load (SWL)	Portable: 286 or 440 lb Fixed: 286, 440, or 600 lb
Noise level	Maximum 54 dBA
Lifting speed	2.2″/sec at 286 lb 1.8″/sec at 440 lb
Battery type	Lithium-ion 25.2V 2500 mAh
Lifts per charge	Approximately 80 cycles of 24" at 200 lb
Battery charging	Full capacity in approximately 2 hours
Life expectancy	10 years with only one battery replacement
Warranty	3-year warranty
Compliance	Meets international certification and compliance standards for medical devices

See our complete line of Silvalea slings at silvaleaslings.com

