

Geo-Mattress® HC

Your Best Bet in Group I Surfaces

Specifically designed to aid in pressure ulcer prevention in the home care setting. Geo-Mattress HC is cut from a single layer of high-quality, combustion modified foam in a fluid-proof, nylon and vinyl cover. Cross-cut top surface is segmented in three distinct zones. Folded and boxed to save freight and storage.

Note: This product does not include Heel Slope® feature.

Item #	Description
SP812-29	Geo-Mattress® HC, 75"L x 35"W x 5"H
SP810-29	Geo-Mattress® HC, 80"L x 35"W x 5"H
SP811-29	Geo-Mattress® HC, 84"L x 35"W x 5"H
Medicare Approved: E0184 (Group I)	
Weight limit: 250 lbs.	Mattress Weight: 17 lbs.
Warranty: 2 years, non-prorated [applies in Home Care use only]	

All models are certified for conformance to NFPA 101 (Life Safety Code) – ASTM E1590; Cal Tech Bulletin #117 and 16 CFR Parts 1632 and 1633 by an independent testing organization.

⚠️ WARNING: Cancer and Reproductive Harm – www.P65Warnings.ca.gov

Description

- ✓ Therapeutic Home Care Foam Mattress
- ✓ Non-powered Pressure Redistribution



Geo-Mattress® 350

High-Value Option Meets Multiple Needs

Segmented top surface, 6" height and fluid-proof cover make the new 350 a great choice for heavier home care patients. Or, an economical yet effective option for long term care residents.


Item #	Description
67535-29	Geo-Mattress 350, 75"L x 35"W x 6"H
68035-29	Geo-Mattress 350, 80"L x 35"W x 6"H
68435-29	Geo-Mattress 350, 84"L x 35"W x 6"H
Medicare Approved: E0184 (Group I)	
Weight limit: 350 lbs.	Mattress Weight: 21 lbs.
Warranty: 5 years, (2 years prorated) [applies in all care settings]	

All models are certified for conformance to NFPA 101 (Life Safety Code) – ASTM E1590; Cal Tech Bulletin #117 and 16 CFR Parts 1632 and 1633 by an independent testing organization.

⚠️ WARNING: Cancer and Reproductive Harm – www.P65Warnings.ca.gov

Description

- ✓ 6" Therapeutic Foam Mattress
- ✓ Non-powered Pressure Redistribution



Gel 80™

Your Best Bet in Group I Surfaces

The Gel-80™, an advanced Group 1 support surface gel/foam overlay, delivers effective pressure redistribution by encapsulating multiple gel bladders within a high-density sculpted support base. Features include a wipe-clean, fluid-impervious fabric top surface, elastic corner straps and one-piece construction. Folded and boxed for easy shipping and storage.

Item #	Description
GEL-80	Gel/foam overlay for hospital-sized twin (80x35) mattress
Medicare Approved: E0185	
Weight limit: 300 lbs.	Dimensions: 78"L x 34"W x 3.5"H
Warranty: 1 year	Overlay Weight: 46 lbs.

Components conform to applicable California Technical Bulletin 117 standards.

Description

- ✓ Gel/Foam Mattress Overlay
- ✓ Non-powered Pressure Redistribution

