

# Seating

Category: Prevention with Positioning

## Geo-Matt® Contour

Exclusive Geo-Matt® anti-shearing layer on top, firmer contoured base below creates outstanding positioning and effective pressure redistribution in one. Base features contoured leg channels, abduction/adduction build-up, and pelvic bucketing to help discourage slouching. Standard cover consists of anti-slip top surface and non-skid bottom. Built-in sleeve in bottom of cover allows use of molded Sling-Fill base (not included-see page 46); 18-month warranty. Height: 3". Weight limit 250 lbs.


**Medicare Code:**  
**Skin Protection—E2603.**

**With Anti-Slip Cover**



**With Deluxe Cover**

rugged Cordura nylon  
with supple, stretch fabric top.



Item #	Description	Qty/Cs
262R	Geo-Matt Contour 16"W x 16"L, durafilm coating, deluxe cover	1
262W	Geo-Matt Contour 18"W x 16"L, durafilm coating, deluxe cover	1
GMC1616-96	Geo-Matt Contour, 16"W x16"L w/ anti-slip cover	6
GMC1616-96-1	Geo-Matt Contour, 16"W x16"L w/ anti-slip cover	1
GMC1816-96	Geo-Matt Contour, 18"W x16"L w/ anti-slip cover	6
GMC1816-96-1	Geo-Matt Contour, 18"W x16"L w/ anti-slip cover	1
GMC1818-96	Geo-Matt Contour, 18"W x18"L w/ anti-slip cover	6
GMC1818-96-1	Geo-Matt Contour, 18"W x18"L w/ anti-slip cover	1
GMC2018-96	Geo-Matt Contour, 20"W x18"L w/ anti-slip cover,	6
GMC2018-96-1	Geo-Matt Contour, 20"W x18"L w/ anti-slip cover	1

Fire barrier versions can be ordered to conform with 2000 NFPA 101 (Life Safety Code) and Cal TB# 133. To order, replace the "-96" suffix with "-29" in item number.

### QUICK FACTS

**Use for:**

- Comfort
- Prevention of breakdown
- Assistance with positioning

**RISK LEVEL**



## Geo-Matt® Wedge Cushion

High resiliency, high density foam wedge has a 4" (front) to 2" (back) taper that helps keep users' hips from sliding forward. Facilitates restraint-free positioning without compromising skin integrity, thanks to its tissue-friendly Geo-Matt® surface geometry. Cover consists of anti-slip top surface and non-skid bottom. Built-in sleeve in bottom of cover allows use of molded Sling-Fill base (not included-see page 46); 18-month warranty. Weight limit 250 lbs.

**Medicare Code:**  
**General Use – E2601.**



Item #	Description	Qty/Cs
GMW1616-96	Geo-Matt Wedge w/ anti-slip cover, 16"W x16"L	6
GMW1616-96-1	Geo-Matt Wedge w/ anti-slip cover, 16"W x16"L	1
GMW1816-96	Geo-Matt Wedge w/ anti-slip cover, 18"W x16"L	6
GMW1816-96-1	Geo-Matt Wedge w/ anti-slip cover, 18"W x16"L	1
GMW1818-96	Geo-Matt Wedge w/ anti-slip cover, 18"W x18"L	6
GMW1818-96-1	Geo-Matt Wedge w/ anti-slip cover, 18"W x18"L	1
GMW2018-96	Geo-Matt Wedge w/ anti-slip cover, 20"W x18"L	6
GMW2018-96-1	Geo-Matt Wedge w/ anti-slip cover, 20"W x 18"L	1

Fire barrier versions can be ordered to conform with 2000 NFPA 101 (Life Safety Code) and Cal TB# 133. To order, replace the "-96" suffix with "-29" in item number.

### QUICK FACTS

**Use for:**

- Comfort
- Prevention of breakdown
- Assistance with positioning

**RISK LEVEL**



**WARNING:** Cancer and Reproductive Harm – www.P65Warnings.ca.gov