


Choice of 6 finishes in 8 profiles




Valley Oak unavailable for Elite profiles  
Bleached Oak also available for Standard profiles

Standard Melamine Profiles



Elite HPL Profiles



# ADVANTAGE

## WEIGH SCALE BED

Span's Weigh Scale low-bed design with advanced clinical positioning and ReadyWide® integrated expansion now available with built-in scale



Position shown:  
Auto Contour


Standard features include:

- Durable, anti-microbial polyester powder-coating of metal components
- Smooth deck design for easier cleaning
- Mobility at any height
- Foot end controller with lockout
- Vascular foot rest with easy ratchet operation
- Beds move vertically with virtually no horizontal motion
- Pendant bed controller
- Advanced mobility as well as floor level locking
- Underbed lighting

Optional features:

- Molded half-head assist rail
- Bed carrier for streamlined bed logistics and maintenance
- Stylish head and foot board options
- Bed extenders

Build your own bed online  
with our interactive  
configurator at  
[www.configurebed.com](http://www.configurebed.com)



# ADVANTAGE

## WEIGH SCALE BED

Height Floor to Mattress Support	8-5/8" – 32-1/4" / 219 mm - 819mm
Deck Angles - Back Elevation	0 – 70°
Deck Angles - Knee Elevation	0 – 30°
Bed Length - Overall Length	82", 86", 90" / 2083 mm, 2184 mm, 2286 mm
Bed Length - Sleep Surface	76", 80", 84" / 1930 mm, 2032 mm, 2134 mm
Bed Width - Sleep Surface	35", 39", 42" / 889 mm, 991 mm, 1067 mm
Weight of Bed (without accessories)	297 lbs / 135 kg
Safe Working Load	500 lbs / 227 kg
Warranty	Mechanical: 5 years Frame: 10 years Welds: 15 years

Item #	Description
MC7WSTML	Weigh Scale Advantage ReadyWide with mobility, floor lock, premium foot end controls and touch screen
MC7WSML	Weigh Scale Advantage ReadyWide with mobility, floor lock & foot end controls
R9938-R	3 position molded pivot assist rail with right bracket
R9938-L	3 position molded pivot assist rail with left bracket
R9939	Set of 3 position molded pivot assist rails
R9308H	Set of molded head end half side rails
R9308H-L	Left side molded head end half side rail
R9308H-R	Right side molded head end half side rail
R9519	Split side rail, 3/4 bed coverage
Q6770	Battery back up - factory installed
Q6740	Floor placement arm
Q6741	Retractable power cord cable reel & placement arm
Q6690	Bed extender to 84"
Q6730	Patient helper, trapeze bar, chain and hand grip
Q6731	Trapeze bed mounting bracket
Q6703	Bed and deck angle indicator
Q6699	Bed transportation system
Q6531	Honeycomb roller bumper kit

Contact your local Span Representative  
or call Beds/Furniture Customer Service:  
**800-268-8671**




Premium Digital  
Weigh Scale Foot  
End Control  
with All-In-One  
Touch Screen  
Operation

# ADVANTAGE

## WEIGH SCALE BED

The Advantage ReadyWide Weigh Scale Bed with advanced features including mobility at any height, easy optional sleep surface expansion to 42" (1067 mm), a multi-function carrier system and a safe working load (SWL) of 500 lbs (227 kg).

Now available with choice of two weigh scale package options:



OR

