Span F450T

Standing Transfer Aid



- Keeps patient actively engaged in the transfer process.
- **Promotes independence,** weight-bearing, and trunk strength.
- High Capacity: 450 lb (205 kg) safe working load.
- Care Flexibility: Underbed clearance allows use with most low beds.







BED















F450T Standing Transfer Aid

The ideal solution for patients with the ability to participate in their own transfers. Makes transport quick and simple while encouraging weightbearing as well as maintenance of trunk and upper body strength.





Swing-away padded seat for secure, comfortable transport



Non-slip footplate provides stable base for transfer



Underbed clearance, allows use with most low beds



Dual, contoured shin pads for added stability

SPECIFICATIONS		
Weight capacity	450 lb / 205 kg	
Overall height	43.2"	
Overall base length	35.5"	
Min. seat height	26"	
Max. seat height	33.5"	
Front-end base height	3.6"	
Internal leg width (open)	32"	
External leg width (closed)	27.7"	
Horizontal handlebar width	20"	
Net Weight	78 lb / 35.2 kg	
Shipping Weight	93 lb / 42 kg	

- Easy base width adjustment accommodates chairs and commodes.
- A convenient toileting alternative vs. use of a wheelchair.
- Locking rear casters add an extra measure of safety.
- Attachment hooks for optional support belt.
- Simple, fast assembly using tools provided.
- Conforms to: ISO 12182:2012.

Item#	Description	
F450T	Sit-To-Stand Patient Lift	
Warranty:	3 years (frame)	

Part of Span's line of Safe Patient Handling and Mobility solutions, including:

F500S Sit-to-Stand Lift

F500P Full Body Patient Lift



Fixed Lithium-ion





Overhead Lifts



Portable Lithium-ion





Floor Lifts