Premium Overbed Table





pictured: flat edge

Span® Premium Overbed Table

Once considered a commodity, overbed tables have become an indispensable component of safe, efficient resident care. This makes their design crucial as well! **Proudly built in North America**, Span tables deliver unmatched, trouble-free value to caregivers and users alike.



pictured: raised edge



Wide, C-shaped base; no center bar to block wheelchair as with typical "H-shaped" bases.



Raised edge helps contain spills, minimizing staff time spent on soiled clothing and wet floors.



Stable, tip-resistant base helps keep table upright; 2" gliding casters.



"Infinite" height adjustment; Spring-assist telescoping column allows smooth raising & lowering.

- Scuff-resistant, anti-microbial polyester powder-coating on base.
- **Seamless top** for superior disinfection (raised edge models only).
- Ships flat: pre-drilled holes for simple, fast assembly.
- Non-prorated warranty: 12 months on all components.
- **Table top dimensions:** 16" Deep x 32" Wide (41cm x 81 cm).
- **Height range:** 28" 45.5" (71 cm 115.5 cm).

Durable Thermoclad DVL

(raised edge models), or low pressure laminate (flat edge) in 6 color options to complement Span's entire line of head/foot boards and case goods.



Choco Cherry



Shiraz Cherry



Summer Flame



Candle Light



Valley Ōak



Hard Rock Maple

Raised Edge	0V11-VLC0	OV11-VLSH	OV11-VLSF	0V11-VLCL		OV11-VLHR
Raised Edge with Vanity	0V12-VLC0	OV12-VLSH	0V12-VLSF	0V12-VLCL		OV12-VLHR
Flat Edge	OC11-MTCOCO	OC11-MTSHSH	OC11-MTSFSF	OC11-MTCLCL	OC11-MTVAVA	OC11-MTHRHR
Flat Edge with Vanity	OC12-MTCOCO	OC12-MTSHSH	OC12-MTSFSF	OC12-MTCLCL	OC12-MTVAVA	OC12-MTHRHR



